

NMR TUBE APPLICATION GUIDE

	<ul style="list-style-type: none"> •Routine •Quick Check •Student Use 	<ul style="list-style-type: none"> •Medium or Heavy Wall •Pressure 	<ul style="list-style-type: none"> •High Resolution •2-D, 3-D, 4D •Long-term Accum. 	<ul style="list-style-type: none"> •Sample Limited •High Filling Factor •Micro Samples 	<ul style="list-style-type: none"> •Low Resolution •Wide Linewidths •Powder Patterns
--	--	--	--	---	---

60 MHz

4mm	405-PS	-	426-PP 427-PP	-	405-PS
5mm	WG-5mm-ECONOMY 512 505-PS	501-PS 502-PS	507-PP 526-PP	529/507-PP 537-PPT	WG-5mm-ECONOMY 512 505-PS
Large Volume	-1PP -3PP	-7PPM -7PPH	-5PP	529A & E / -3PP -7PPT or -7PPXT	-1PP

80 - 100 MHz

4mm	-	-	-	-	-
5mm	506-PP 507-PP	502-PP 504-PP	526-PP 527-PP 528-PP	529-A/507-PP 537-PPT	505-PS 506-PP
Large Volume	-1PP -3PP	-7PPM -7PPH	-5PP -7PP	529-A & E / -5PP -7PPT or -7PPXT	-3PP

150 - 250 MHz

4mm	-	-	-	-	-
5mm	507-PP 526-PP	522-PP 524-PP	527-PP 528-PP 535-PP	529-A/528-PP 540-PPT	506-PP 507-PP
Large Volume	-3PP -5PP	-7PPM -7PPH	-5PP -7PP	529A & E / -5PP -7PPT or -7PPXT	-5PP

270 - 400 MHz

4mm	-	-	-	-	-
5mm	526-PP 527-PP	522-PP 524-PP	528-PP 535-PP	529-A/528-PP 540-PPT	506-PP 507-PP
Large Volume	-5PP -7PP	-7PPM -7PPH	-7PP	529A & E / -5PP -7PPT or -7PPXT	-5PP

400 - 500 MHz

4mm	-	-	-	-	-
5mm	526-PP 527-PP	522-PP 524-PP	528-PP 535-PP	529-A / 528-PP 540-PPT	506-PP 507-PP
Large Volume	-5PP -7PP	-7PPM -7PPH	-7PP	529A & E / -5PP -7PPT or -7PPXT	-5PP

600 - 750 MHz +

4mm	-	-	-	-	-
5mm	526-PP 527-PP	522-PP 524-PP	528-PP 535-PP	529-A/528-PP 540-PPT	506-PP 507-PP
Large Volume	-	-	-	-	-

弊社在庫品